



Acropora DNA extraction with Qiagen DNeasy tissue kit V.2

Iliana Baums¹, Sheila A Kitchen¹

¹Pennsylvania State University

Version 2

May 18, 2020

1

Works for me

dx.doi.org/10.17504/protocols.io.bgjqjumw

Astrangia

Coral Genotyping Protocols

Iliana Baums

ABSTRACT

DNA extraction protocol for Acropora or other coral tissue based on Qiagen DNeasy kit. This extraction protocol works well for the Acropora SNPchip and other coral genotyping applications (such as microsatellite genotyping). It preferentially extracts coral host DNA but some *Symbiodiniaceae* DNA will be present.

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Baums IB, Hughes CR, Hellberg MH (2005) Mendelian microsatellite loci for the Caribbean coral *Acropora palmata*. Mar Ecol Prog Ser 288:115-127. doi:10.3354/meps288115.

ATTACHMENTS

[Acroporid tissue DNA extraction.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.bgjqjumw

PROTOCOL CITATION

Iliana Baums, Sheila A Kitchen 2020. Acropora DNA extraction with Qiagen DNeasy tissue kit. **protocols.io** dx.doi.org/10.17504/protocols.io.bgjqjumw

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Baums IB, Hughes CR, Hellberg MH (2005) Mendelian microsatellite loci for the Caribbean coral *Acropora palmata*. Mar Ecol Prog Ser 288:115-127. doi:10.3354/meps288115.

KEYWORDS

Acropora, DNA extraction, SNPchip, STAGdb

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 18, 2020

LAST MODIFIED

May 26, 2020

PROTOCOL INTEGER ID

37200

GUIDELINES

This protocol works well on coral tissue preserved in 95% non-denatured Ethanol. Do not overload the columns. A few polyps is enough. Make sure your centrifuge can spin high enough to complete this protocol before you start. Both the 96-well plate version and the single 1.5 ml tube version of the QIAGEN kit work. For the 1.5ml tube version, we like the electronic multi-channel MATRIX pipettor with flexible spacing to go from 1.5 ml tubes to PCR tubes/plates. Avoid defrost cycles of samples and extracted DNA.

MATERIALS

NAME	CATALOG #	VENDOR
QIAgen DNeasy Blood and Tissue Kit, 50 rxn	69504	Qiagen
Proteinase K, 2mL	19131	Qiagen
Incubator		
Razor blade		
Centrifuge	5415 D	Eppendorf Centrifuge

Sample preparation

- 1 Use sterile technique to cut 3-4 polyps (skeleton with tissue) off the coral fragment with a razor blade.
- 2 Transfer the pieces to a 1.5 ml tube/96-well plate.
- 3 Add 180 ul of Buffer ATL and 20 ul of proteinase K (from Qiagen kit) to each tube/well.
- 4 If using the tubes, invert the tubes several times to gently mix the reagents. *Do not vortex*
- 5 Incubate the mixture overnight at 56 °C.
- 6 a. The mixture will have a slight yellow hue after lysis from the algal symbionts.
- 7 b. The tubes or plate can be placed on a rocking platform for gentle agitation.

Qiagen DNeasy Blood and Tissue Kit protocol

- 8 Follow Qiagen's protocol with the following changes:
Following step 6 of the spin-column protocol, we perform the optional centrifugation for 1 min at 20,000 xg to dry the column membrane.

Elution 2x

- 9 Elution 1- Add 100 ul of low TE buffer directly to the column membrane, incubate for 1 min and centrifuge as recommended.
a. This elution typically contains smaller fragmented DNA.

- 10 Elution 2 – Move the column to a new 1.5 ml tube. Repeat elution as above, but increase incubation of 100 ul of low TE to 5 min.
a. This elution contains high molecular weight DNA and is used for downstream analyses.
b. Store at -20 °C.

Check extraction quality

- 11 Run DNA on a 1.5% agarose gel to check extraction quality of elution 1 and 2. Look for one high molecular weight band at the top of the gel with little smearing.