Pham, H., Howe, J., Ghosh, A., Melton, M., Kunwar, D., Peterson, E., & Datye, A. (2018). Using a Combination of HAADF and SE Imaging to Locate Pt Nanoparticles within a Mesoporous Silica Diesel Oxidation Catalyst. *Microscopy and Microanalysis*, 24(S1), 1700-1701. doi:10.1017/S143192761800898X

Using a Combination of HAADF and SE Imaging to Locate Pt Nanoparticles within a Mesoporous Silica Diesel Oxidation Catalyst

Hien N. Pham¹, Jane Howe², Arnab Ghosh¹, Matthew Melton¹, Deepak Kunwar¹, Sivakumar R. Challa¹, Eric J. Peterson¹, and Abhaya K. Datye¹.

The USDRIVE Low Temperature Oxidation Catalyst Test Protocol recommends an accelerated aging test for Pt diesel oxidation catalysts (DOCs) to be performed at 800 °C in flowing air. However, at such high temperatures, Pt sinters readily to form large particles due to significant vapor phase transport of PtO₂. One approach we used to lower the vapor pressure of PtO₂ was by alloying Pt with Pd, which helped to reduce the rate of Pt sintering [1]. A second approach was to trap the mobile Pt species with an oxide such as PdO that reacts to form Pt-Pd bimetallic particles [2], or using CeO₂, which traps Pt in atomically dispersed form [3]. Here we explore an alternative strategies, such as the use of pore confinement, to determine if it is possible to slow the vapor phase transport of PtO₂. Since PtO₂ vapor has negligible binding to silica [4], it provides the ideal support to study the role of pore structure. In future work, we will explore different pore structures and connectivity; for this work we focused on the approach for imaging nanoparticles in mesoporous SBA-15 silica.

The Pt diesel oxidation catalysts were synthesized via incipient wetness impregnation using mesoporous silica SBA-15 as the support, and tetra amine platinum (II) nitrate as the precursor (nominal loading 1.5 wt% Pt). The Pt DOCs were loaded inside a tube furnace and aged in air (50 cc/min) at 800 °C for up to 5 h. HAADF-STEM was performed using a JEOL 2010F at 200 kV, equipped with an Oxford Aztec energy dispersive X-ray spectroscopy (EDS) system. To determine whether some of the Pt has remained inside the pores of SBA-15 after aging, simultaneous dark field (DF) and secondary electron (SE) modes were performed using a Hitachi probe aberration-corrected HF-5000 operated at 60 kV.

Figure 1(a) shows the pore size distribution SBA-15 silica. We see a narrow pore size distribution, consistent with a well-ordered structure, and an average pore size of 6 nm. After depositing Pt, the asprepared Pt/SBA-15 (Figure 1(b)) catalyst contains highly dispersed Pt particles ~1-2 nm in size. Aging in air at 800 °C for 5 h causes severe sintering of the catalyst, causing a bimodal distribution with the coexistence of small (~2 nm) and large (10 nm - 100 nm) Pt particles in Pt/SBA-15 (Figure 1(b)). After aging, the DF image shows only the outlines of the large, faceted Pt particles whereas the SE images show the surface morphology, indicating that the large Pt particles are present on the outer surface of SBA-15 (Figure 1(c)). The DF image also shows the presence of small Pt particles but they are not observed in the SE image, indicating that the small Pt particles are retained inside the pores of SBA-15 (Figure 1(c)). By performing EDS from regions that only contain small particles, we concluded that for Pt/SBA-15, 0.5 wt% Pt is retained in the form of small particles. Figure 2 shows the DF and SE images of Pt/SBA-15 aged at 800 °C for (a) 1, (b) 3, and (c) 5 h. The pores of SBA-15 help to confine and retain the Pt with the result that only few large Pt particles are seen after 1 h and aging. The number of large Pt particles progressively increases as the catalyst is aged, and after 5 h and significant fraction of the Pt is now located on the external surface of the silica support.

^{1.} Chemical & Biological Engineering and Center for Micro-engineered Materials, University of New Mexico, Albuquerque, NM 87110, USA


² Hitachi High-Technologies America Inc., Clarksburg, MD 20871, USA

Pham, H., Howe, J., Ghosh, A., Melton, M., Kunwar, D., Peterson, E., & Datye, A. (2018). Using a Combination of HAADF and SE Imaging to Locate Pt Nanoparticles within a Mesoporous Silica Diesel Oxidation Catalyst. *Microscopy and Microanalysis*, *24*(S1), 1700-1701. doi:10.1017/S143192761800898X

The significance of this work is that it presents a strategy to control the rate of Pt sintering. The cylindrical pores of SBA-15 help to retain some of the Pt within the pores, which then form small particles that co-exist with the larger particles that form outside the pores. The SE images are critical to understanding the location of the Pt and the role of pore structure in improving the sinter resistance [5].


References:

- [1] C. Carrillo et al, Applied Catalysis B: Environmental 218 (2017), p. 581.
- [2] C. Carrillo et al, Journal of Physical Chemistry Letters 5 (2014), p. 2089.
- [3] J. Jones *et al*, Science **353** (2016), p. 150.
- [4] P.N. Plessow and F. Abild-Pedersen, ACS Catalysis 6 (2016), p. 7098.



[5] The authors acknowledge funding from the NSF GOALI Grant CBET-1438765.

Figure 1. (a) Pore size distribution of SBA-15 support. HAADF-STEM (200 kV) images of (b) Pt/SBA-15, as prepared and after aging in air at 800 °C for 5 h. DF and SE images (60 kV) of (c) Pt/SBA-15 after aging for 5 h.



Pham, H., Howe, J., Ghosh, A., Melton, M., Kunwar, D., Peterson, E., & Datye, A. (2018). Using a Combination of HAADF and SE Imaging to Locate Pt Nanoparticles within a Mesoporous Silica Diesel Oxidation Catalyst. *Microscopy and Microanalysis*, *24*(S1), 1700-1701. doi:10.1017/S143192761800898X

Figure 2. Dark field (DF) and secondary electron (SE) images (60 kV) of Pt/SBA-15 aged in air 800 °C for (a) 1, (b) 3, and (c) 5 h. The Pt initially in the pores (see Figure 1(b)) progressively comes out of the pores, forming large Pt particles located on the surface of the silica support.