



# Biochemical Approaches for Understanding Iron–Sulfur Cluster Regeneration in *Escherichia coli* Lipoyl Synthase During Catalysis

Erin L. McCarthy\*, Squire J. Booker\*,†,‡,¹

\*Department of Biochemistry and Molecular Biology, The Pennsylvania State University, University Park, PA, United States

†Department of Chemistry, The Pennsylvania State University, University Park, PA, United States

‡The Howard Hughes Medical Institute, The Pennsylvania State University, University Park, PA, United States

¹Corresponding author: e-mail address: [squire@psu.edu](mailto:squire@psu.edu)

## Contents

1. Introduction	218
2. Overproduction and Isolation of <i>E. coli</i> LipA and <i>E. coli</i> NfuA	220
2.1 Gene Cloning and Expression Strategy	220
2.2 Protocol for Expression of the <i>lipA</i> and <i>nfuA</i> Genes	221
2.3 Purification of <i>E. coli</i> LipA and NfuA	222
3. Initial Characterization of As-Isolated Protein	224
3.1 Protocol for Obtaining Amino Acid Correction Factor	226
3.2 Colorimetric Iron Analysis Protocol	226
3.3 Colorimetric Sulfide Analysis Protocol	227
4. Methods for Isotopically Labeling Iron–Sulfur Clusters to Monitor Cluster Transfer	228
4.1 Protocol for the Expression of the <i>nfuA</i> Gene in Apo Form	228
4.2 Synthesis of <sup>34</sup> S-Labeled Sulfide From <sup>34</sup> C-Labeled Elemental Sulfur	229
4.3 Reconstitution of Apo NfuA Using Na <sub>2</sub> <sup>34</sup> S	230
5. Strategies for Analyzing Protein Complex Formation	230
5.1 Protocol for Performing Analytical Molecular-Sieve Chromatography	232
6. Monitoring Iron–sulfur Cluster Transfer Using Quantitative Analysis	233
6.1 Sample Preparation for Analysis of Reaction Products by LC-MS to Assess the Effect of NfuA on LipA Activity	234
6.2 Sample Preparation for Analysis of Reaction Products in the Presence of a Competitive Chelator or Exogenous <sup>34</sup> C-Labeled Sulfide to Ensure Direct Cluster Transfer From NfuA to LipA	234
6.3 Protocol for the Quantification of Reaction Products by LC-MS	235
7. Conclusions	237
Acknowledgments	237
References	237
Further Reading	239

## Abstract

Lipoyl synthase (LipA in bacteria) is a radical S-adenosylmethionine (SAM) enzyme that catalyzes the second step of the de novo biosynthesis of the lipoyl cofactor: the insertion of sulfur at C6 and C8 of a pendant octanoyl chain. In addition to the [4Fe—4S] cluster that is characteristic of the radical SAM (RS) enzymes, LipA contains a second [4Fe—4S] cluster that, though controversial, has been proposed to be degraded during turnover to supply the inserted sulfur atoms. A consequence of this proposed role is that the destruction of its iron–sulfur cluster renders the enzyme in an inactive state. Recently, it was shown that *Escherichia coli* proteins NfuA or IscU can confer catalytic properties to *E. coli* LipA in vitro. In this chapter, we present methods for characterizing LipA and analyzing its activity in vitro, and provide strategies to monitor the pathway for the regeneration of LipA's auxiliary cluster by *E. coli* iron–sulfur carrier protein NfuA.



## 1. INTRODUCTION

Lipoic acid is an eight carbon straight-chain fatty acid containing sulfur atoms at C6 and C8 (Reed, Busk, Gunsalus, & Schnakenberg, 1951; Reed et al., 1953). It plays a central role in multienzyme complexes involved in the oxidative decarboxylation of  $\alpha$ -keto acids, the breakdown of glycine, and the degradation of acetoin in certain types of bacteria (Reed, 1974). In this capacity, it functions as a redox-active cofactor due to its ability to convert between states containing free thiol groups and a dithiolane ring and as a shuttle to move reaction intermediates between successive active sites in multienzyme complexes. In its role as a cofactor, lipoic acid is covalently attached through an amide linkage to a specific lysyl residue of a lipoyl carrier protein (LCP), one of the subunits of the multienzyme complexes (Cronan, 2016).

Early studies identified the pathway for the de novo biosynthesis of the lipoyl cofactor in *Escherichia coli* and outlined the roles of two dedicated enzymes: an octanoyltransferase, LipB, and a radical S-adenosylmethionine (SAM) enzyme, LipA (Cronan, 2016; Miller et al., 2000; Zhao, Miller, Jiang, Marletta, & Cronan, 2003). In an offshoot from fatty acid biosynthesis, LipB transfers an octanoyl group from an acyl carrier protein to an LCP, which generates the substrate for LipA. LipA subsequently inserts sulfur atoms at C6 and C8 to generate the lipoyl cofactor. The inability to synthesize the lipoyl cofactor is embryonic lethal in mice, and defects in its biosynthesis result in disastrous biochemical outcomes in humans, such as seizures, developmental delays, and death (Landgraf, McCarthy, & Booker, 2016; Yi & Maeda, 2005). Further, lipoic acid supplementation has been shown to be therapeutic for a wide range of clinical ailments

associated with oxidative stress, including diabetes and Alzheimer's disease (Borcea et al., 1999; Holmquist et al., 2007).

There are several challenges associated with the biosynthesis of the lipoyl cofactor, despite its structural simplicity. The first obstacle is removal of hydrogens from C6 and C8 of the octanoyl chain, which LipA overcomes through the use of radical chemistry. As a radical SAM (RS) enzyme, LipA coordinates a  $[4\text{Fe}-4\text{S}]^{2+}$  cluster that is absolutely required for its function. In its reduced  $[4\text{Fe}-4\text{S}]^+$  form, the iron-sulfur (Fe/S) cluster can inject an electron into the sulfonium of SAM, inducing its cleavage to L-methionine, and a 5'-deoxyadenosyl 5'-radical (5'-dA $\cdot$ ). This high energy radical species is a potent oxidant capable of cleaving C—H bonds displaying homolytic bond-dissociation energies of 98 kcal/mol or greater.

The second outstanding question is the source of the sulfur used for the sulfur insertion. In addition to the  $[4\text{Fe}-4\text{S}]$  cluster common to all RS enzymes, LipA coordinates a second “auxiliary” cluster that appears to be the source of the inserted sulfur atoms (Lanz & Booker, 2015). As such, the  $[4\text{Fe}-4\text{S}]$  auxiliary cluster is destroyed during catalysis, which is supported by an overwhelming amount of biochemical evidence (Cicchillo & Booker, 2005; Cicchillo, Iwig, et al., 2004; Cicchillo, Lee, et al., 2004; Douglas, Kriek, Bryant, & Roach, 2006; Harmer et al., 2014; Lanz et al., 2014; McLaughlin et al., 2016). For example, LipA can only undergo up to one turnover in vitro in the absence of exogenously added sulfur. Further, when LipA is isolated from *E. coli* cultures using  $^{34}\text{S}$  as the sole sulfur source and then subsequently used in in vitro assays, the lipoyl product contains almost exclusively  $^{34}\text{S}$  at C6 and C8 (Cicchillo & Booker, 2005).

Due to the sacrificial role of its auxiliary cluster, LipA typically catalyzes less than one equivalent of product formation in in vitro assays. Recently, the *E. coli*  $[4\text{Fe}-4\text{S}]$  cluster carrier protein NfuA was shown to regenerate LipA's auxiliary cluster to support catalysis in vitro (McCarthy & Booker, 2017). Remarkably, it accomplishes this feat at a rate that does not limit the rate of catalysis. In the absence of NfuA, the scaffold protein IscU can fulfill a similar role. Herein, we describe methods for isolating and characterizing LipA and monitoring the regeneration of its auxiliary cluster in vitro. This chapter aims to complement other previously published *Methods in Enzymology* papers that highlight the characterization of RS enzymes that contain multiple  $[\text{Fe}-\text{S}]$  clusters (Lanz et al., 2012) and RS enzymes that contain cobalamin cofactors (Błaszczuk, Wang, & Booker, 2017), with additional insight into analysis of Fe—S cluster regeneration in vitro in RS enzymes that consume their clusters during catalysis.



## 2. OVERPRODUCTION AND ISOLATION OF *E. COLI* LipA AND *E. COLI* NfuA

The generation of LipA protein of high purity and activity has been notoriously difficult due to issues of solubility and the air sensitive nature of its two [4Fe–4S] clusters. Initial characterizations resulted in protein that had poor activity most likely due to incomplete cofactor incorporation which is essential for the enzyme's activity. With regard to NfuA, the *E. coli* protein had been typically isolated aerobically and then subjected to a chemical reconstitution procedure to insert its iron–sulfur (Fe/S) cluster, leaving it unclear as to what the *in vivo* cluster configuration is (Angellini et al., 2008; Boutigny et al., 2013; Py et al., 2012). Following the methods detailed herein, we have generated robust amounts of LipA and NfuA proteins that are >90% homogenous (Figs. 1 and 2). Further, the proteins contain intact [4Fe–4S] clusters in its As-isolated state without the need for chemical reconstitution allowing for extensive biochemical analysis

### 2.1 Gene Cloning and Expression Strategy


The genes encoding LipA and NfuA were amplified from *E. coli* genomic DNA (strain W3110) using the polymerase chain reaction (PCR) and inserted into the pET28a vector using the *NdeI*/*EcoRI* restriction sites and *HindIII*/*XhoI* restriction sites, respectively (Cicchillo, Iwig, et al., 2004; Cicchillo, Lee, et al., 2004; McCarthy & Booker, 2017). The resulting constructs allowed for overexpression of the *lipA* and *nfuA* genes under the control of the strong T7 promoter. Further, these constructs conferred kanamycin resistance for selection *in vivo* and encoded an N-terminal hexahistidine (His<sub>6</sub>) tag, which provided a simple purification strategy. A second plasmid, pDB1282, encodes the *isc* operon from *Azotobacter vinelandii* and confers ampicillin resistance, making it compatible for coexpression with the pET28a-*lipA* and pET28a-*nfuA* constructs. *E. coli* BL21(DE3) cells, which are commonly used with the pET systems, were cotransformed with either the pET28a-*lipA* or pET28a-*nfuA* construct and pDB1282 for gene expression.

*E. coli* LipA and *E. coli* NfuA are overproduced in *E. coli* cultured in M9 minimal medium supplemented with glucose as a carbon source and ammonium chloride as a nitrogen source, and cultures are grown with limited aeration. During the overproduction, iron is not typically added to the cultures until gene expression is induced at  $\sim\text{OD}_{600}=0.6$ . We have found that in some cases, iron starvation can induce cells to release siderophores,


which have an extremely high affinity for iron and can increase iron–sulfur cluster incorporation in some proteins, such as LipA.

## 2.2 Protocol for Expression of the *lipA* and *nfuA* Genes

1. Cotransform *E. coli* BL21(DE3) cells with pDB1282 and either pET28a-*lipA* or pET28a-*nfuA* using standard protocols. Use the antibiotic markers within the plasmids to select for colonies containing both pDB1282 and pET28a-*lipA* or pET28a-*nfuA* constructs by plating transformants on Luria-Bertani (LB) plates containing kanamycin (50 µg/mL) and ampicillin (100 µg/mL).
2. Select a single colony and inoculate 200 mL LB media containing kanamycin (50 µg/mL) and ampicillin (100 µg/mL). Incubate at 37°C with shaking at 250 rpm for ~14 h. Simultaneously, place four 6 L flasks containing 4 L M9 minimal media in a shaker overnight at 37°C in order to warm culture media.
3. Immediately prior to inoculation of large flasks, add kanamycin (50 µg/mL) and ampicillin (100 µg/mL).
4. Inoculate each 6 L flask with 10 mL of the 200 mL overnight LB culture and incubate the cultures at 37°C with shaking at 180 rpm.
5. Once the cells reach an  $OD_{600}=0.3$  (~5 h), add 0.2% (w/v) L-arabinose to induce expression of genes on pDB1282. Take a 1 mL sample of each flask directly before arabinose induction (preinduction sample) for SDS-PAGE analysis.
6. At an  $OD_{600}=0.6$ , place the flasks in an ice bath for ~1 h and add  $FeCl_3$  to a final concentration of 50 µM. Adjust the temperature of the incubator to 18°C. Take a sample of each flask (postarabinose sample) for SDS-PAGE analysis, normalizing the number of cells to the preinduction sample.
7. Induce gene expression by adding 200 µM isopropyl-β-D-thiogalactopyranoside (IPTG) and allow expression to proceed for ~18 h at 18°C with shaking at 180 rpm.
8. The following day, take a sample of each flask (post-IPTG sample) for SDS-PAGE analysis, normalizing the number of cells to the postarabinose sample. Harvest the bacterial cells by centrifugation at  $7500 \times g$  for 12 min. Flash-freeze the resulting cellular pellet (~3 g/L) in liquid nitrogen and store in liquid nitrogen until further use.
9. Analyze the overproduction of the desired protein by SDS-PAGE (12.5%) using the samples collected during the growth (preinduction, postarabinose induction, post-IPTG induction).



**Fig. 1** SDS-PAGE analysis of (A) *E. coli* LipA overproduction and (B) *E. coli* NfuA overproduction. The lanes are as follows: (1) Molecular Weight Ladder, (2) Uninduced cells, (3) Arabinose-induced cells, (4) IPTG-induced cells. The *arrows* point to the protein of interest.



**Fig. 2** SDS-PAGE analysis of (A) *E. coli* LipA isolation and (B) *E. coli* NfuA isolation. The lanes are as follows: (1) Molecular Weight Ladder, (2) Pellet, (3) Supernatant, (4) Flow-through, (5) Wash 1, (6) Wash 2, (7) Pure protein. The *arrows* point to the protein of interest.

### 2.3 Purification of *E. coli* LipA and NfuA

Isolating proteins that contain Fe/S clusters are challenging due to their oxygen-sensitive nature. Oxygen contamination can lead to the degradation of the clusters, resulting in proteins with poor stability and enzymatic activity. As such, careful handling of the protein during all purification steps is vital to ensure that the quality of the protein does not suffer because of oxidative damage. Ideally, the isolation and characterization of Fe/S cluster-containing proteins should be performed in an anaerobic chamber.

The Fe/S clusters coordinated by LipA and NfuA are especially sensitive to oxygen, so maintaining an anaerobic environment is essential. LipA is also sensitive to repeated freeze–thaw cycles, so the purification should be carried out to completion without interruption using chilled buffers, and the protein should be stored in small aliquots. The general procedure for isolating Fe—S cluster-containing proteins has been previously described in detail (Lanz et al., 2012).

1. Autoclave plastic centrifuge tubes and cycle into anaerobic chamber while hot. Allow the tubes to equilibrate in the chamber several days before purification.
2. In the anaerobic chamber, prepare 250 mL of each of the buffers needed for the purification using oxygen-free water: lysis buffer (50 mM HEPES, pH 7.5, 300 mM KCl, 20 mM imidazole, 10 mM 2-mercaptoethanol), wash buffer (50 mM HEPES, pH 7.5, 300 mM KCl, 45 mM imidazole, 10 mM 2-mercaptoethanol, 10% glycerol), elution buffer (50 mM HEPES, pH 7.5, 300 mM KCl, 500 mM imidazole, 10 mM 2-mercaptoethanol, 10% glycerol), and storage buffer (50 mM HEPES, pH 7.5, 300 mM KCl, 5 mM dithiothreitol, 20% glycerol).
3. Resuspend the cell paste (~60 g) in 150 mL lysis buffer supplemented with lysozyme (1 mg/mL) and DNase I (0.1 mg/mL) in a metal sonication cup. Once the solution is homogeneous, allow it to stir for 30 min at room temperature and then add one EDTA-free Sigmafast protease inhibitor tablet.
4. Place the sonication cup in an ice-water bath on a stir-plate and slowly stir.
5. Lyse the cells by sonic disruption using a Fisher 550 Sonic Dismembrator by sonicating for 1 min at 30% output (setting of 7) followed by resting for 8 min to allow for temperature reequilibration. Repeat for a total of five cycles. The cell lysate should become noticeably darker as the cells burst, releasing Fe—S containing proteins into solution.
6. Following sonication, pour the cell lysate into the plastic centrifuge tubes that were previously autoclaved and allowed to equilibrate in the anaerobic chamber. Seal the juncture between the cap and tube tightly with vinyl tape to avoid oxygen contamination.
7. Remove the tightly sealed tubes from the glove box and centrifuge at  $45,000 \times g$  for 1 h at 4°C.
8. Cycle the centrifuge tubes back into the anaerobic glove box and pour the supernatant into a 250 mL flask. Take a 50  $\mu$ L sample of the

supernatant and scrape a small amount of the pellet into an Eppendorf tube for SDS-PAGE analysis.

9. Equilibrate a column containing 50 mL Ni-NTA resin with 100 mL chilled lysis buffer.
10. After equilibration, load the supernatant onto the column. Collect the flow-through in a separate container and take a 50  $\mu$ L sample for SDS-PAGE analysis.
11. Wash the column with 100 mL chilled wash buffer and take a 50  $\mu$ L sample for SDS-PAGE analysis. Repeat the wash step and collect another 50  $\mu$ L sample for SDS-PAGE analysis.
12. Elute the desired protein using chilled elution buffer. Allow  $\sim$ 15 mL of elution buffer to flow over the column before collecting the protein, which is when the dark brown colored band will begin to elute. Stop collecting once the color of the eluate returns to clear.
13. Concentrate the protein in an Amicon stir-cell fitted with an YM-10 membrane (10 kDa molecular weight cutoff) in the anaerobic glove box. It is vital to keep the stirring apparatus on ice throughout the entire concentration step.
14. Once the protein has been concentrated to  $\sim$ 2.5 mL, exchange it into chilled storage buffer using a PD-10 column. Collect 3.5 mL and aliquot into Eppendorf tubes.
15. Take a 5  $\mu$ L aliquot of the final purified protein for SDS-PAGE analysis.
16. Remove the protein aliquots from the anaerobic chamber and immediately flash-freeze in liquid nitrogen. Store the protein aliquots in liquid nitrogen until use, avoiding freezing and thawing, which leads to activity loss.
17. Analyze the purification by SDS-PAGE (12.5%) using the samples collected of the cellular pellet, supernatant, flow-through, wash, and final protein.



### **3. INITIAL CHARACTERIZATION OF AS-ISOLATED PROTEIN**

Following the isolation, routine characterizations should be performed to determine the yield and to obtain a qualitative estimate of the cofactor incorporation in the protein. These routine analyses provide a straightforward method to monitor consistency among enzyme preparations. UV-visible spectroscopy can be used to analyze [4Fe-4S] cluster-containing proteins which have distinct features at 330 and 400 nm (Fig. 3). Protein






Fig. 3 UV-visible spectrum of (A) As-isolated *E. coli* LipA and (B) As-isolated *E. coli* NfuA.




Fig. 4 (A) Iron analysis standard curve. (B) Sulfide analysis standard curve.

concentration is determined via the Bradford method (Bradford, 1976) using bovine serum albumin (fraction V) as a standard. However, because the Bradford method is biased by amino acid composition, a correction factor should be obtained through amino acid analysis to account for differences in amino acid sequence. This step was especially important in the analysis of LipA, which has a correction factor of nearly 1.5. In this case, its amino acid composition significantly changes the protein concentration empirically determined by the Bradford method, which has further implications in determining its activity and the stoichiometry of its Fe—S cluster cofactors. In addition to protein concentration, colorimetric iron and sulfide analyses should be performed (Beinert, 1978, 1983) (Fig. 4). It is important to note

that these methods simply estimate the amount of iron and sulfide in solution; they do not provide evidence for the biologically relevant form of the cofactor or its coordination environment. For a more advanced understanding, Mössbauer and EPR spectroscopies should be used, as detailed in other methods papers (Lanz et al., 2012; Pandelia, Lanz, Booker, & Krebs, 2015).

### 3.1 Protocol for Obtaining Amino Acid Correction Factor

1. Contact a facility capable of performing amino acid analysis, such as the University of California-Davis proteomics facility.
2. Ensure that the protein you wish to analyze is >95% homogenous.
3. Using a NICK column (GE healthcare) exchange 100  $\mu$ L of the protein ( $\sim$ 40 mg/mL) into 10 mM NaOH in order to remove the Fe—S cluster.
4. Perform three Bradford assays each with an independent standard curve to determine the protein concentration.
5. Average the three values together and record as the concentration determined via the Bradford method.
6. Carefully place 100  $\mu$ L of the protein in three Corning Pyrex test tubes (10  $\times$  75 mM) and seal the top with parafilm.
7. Use a needle to add a hole to the top of each of the glass tubes and lyophilize the protein samples in vacuo.
8. Ship to the proteomics facility in order to determine the true concentration of the sample based on its amino acid composition.
9. Average the three values of the concentration determined by amino acid analysis.
10. Generate the amino acid correction factor by dividing the concentration determined by amino acid analysis by the concentration determined by the Bradford method.
11. Apply this ratio in every subsequent Bradford analysis performed on the protein to ensure an exact determination of the concentration.

### 3.2 Colorimetric Iron Analysis Protocol

1. Prepare the following solutions: Reagent A (1.35 g SDS, 450  $\mu$ L saturated iron-free ammonium acetate, in a final volume of 30 mL dH<sub>2</sub>O); Reagent B (270 mg ascorbic acid, 9 mg sodium metabisulfite, 450  $\mu$ L saturated iron-free ammonium acetate, in a final volume of 6 mL dH<sub>2</sub>O); and Reagent C (18 mg ferene in a final volume of 1 mL dH<sub>2</sub>O).

2. In the anaerobic chamber, dilute the protein to 10  $\mu\text{M}$  in 50 mM NaOH.
3. Aliquot into four 300  $\mu\text{L}$  samples and remove the samples from the anaerobic chamber.
4. Prepare a 100  $\mu\text{M}$  iron standard stock using Claritas PPT iron standard (Fisher Scientific).
5. Dilute the 100  $\mu\text{M}$  iron standard to the following concentrations in a final volume of 300  $\mu\text{L}$ : 0, 1.66, 3.33, 16.6, 33.3, 50, 66.6, 83.3, and 100  $\mu\text{M}$ .
6. Add 300  $\mu\text{L}$  of Reagent A to each sample and standard and mix.
7. Add 300  $\mu\text{L}$  Reagent B to each sample and standard and mix.
8. Incubate the samples and standards in a 37°C water bath for 15 min.
9. Add 15  $\mu\text{L}$  of Reagent C to each sample and standard. Mix and then centrifuge at  $14,000 \times g$  for 5 min.
10. Measure the absorbances of the standards at 592 nm using the first standard (0  $\mu\text{M}$ ) as a blank.
11. Measure the absorbances of each of the four samples.
12. Plot the concentration vs absorbance of each of the standards and use the linear equation to calculate the concentrations of each of the four samples (Fig. 4).
13. Average the concentration of iron in each of the four samples and normalize to the protein concentration in order to determine the amount of iron per protein.

### 3.3 Colorimetric Sulfide Analysis Protocol

1. Prepare the following solutions anaerobically: 1% zinc acetate (w/v%); 3 M NaOH; 0.2% dimethyl-4-phenylenediamine (DMPD) in 5 N HCl; 23 mM  $\text{FeCl}_3$  in 1.2 N HCl; and 200  $\mu\text{M}$   $\text{Na}_2\text{S}$ .
2. Complete steps 3–10 anaerobically.
3. Dilute the protein to 10  $\mu\text{M}$  in 50 mM NaOH (prepared in deoxygenated water) to a final volume of 1000  $\mu\text{L}$ .
4. Aliquot 100  $\mu\text{L}$  of the 10  $\mu\text{M}$  protein solution into five Eppendorf tubes.
5. Create a standard curve using the 200  $\mu\text{M}$   $\text{Na}_2\text{S}$  solution by making 100  $\mu\text{L}$  standards at the following concentrations: 0, 10, 20, 30, 40, 50, 60, and 80  $\mu\text{M}$ .
6. Add 300  $\mu\text{L}$  of the 1% zinc acetate solution to each of the standards and samples.
7. Add 15  $\mu\text{L}$  of 3 M NaOH to each of the standards and samples. Mix by inversion.

8. Let stand for 2 h to allow zinc sulfide to form and settle to the bottom of the tubes.
9. Underlay the standards and samples with 75  $\mu$ L of the 0.2% DMPD solution.
10. Quickly add 2  $\mu$ L of 23 mM FeCl<sub>3</sub> to each sample and mix by inversion.
11. Remove standards and samples from the anaerobic chamber and centrifuge at  $14,000 \times g$  for 15 min.
12. Remove 400  $\mu$ L of each standard and sample and place in unused Eppendorf tubes.
13. Add 600  $\mu$ L 2 mM EPPS, pH 8.0, to each standard and sample.
14. Let stand for 30 min at room temperature and then measure the absorbance at 666 nm.
15. Plot the concentration vs the absorbance from each standard. Use the linear equation to determine the concentration of sulfide in each sample and normalize the sulfide to the protein concentration to determine the amount of sulfide per protein (Fig. 4).



## **4. METHODS FOR ISOTOPICALLY LABELING IRON–SULFUR CLUSTERS TO MONITOR CLUSTER TRANSFER**

### **4.1 Protocol for the Expression of the *nfuA* Gene in Apo Form**

1. Introduce the pET28a-*nfuA* construct into BL21(DE3) cells using standard transformation procedures. Spread cells on an LB plate supplemented with 50  $\mu$ g/mL kanamycin.
2. The following day, use a single colony to inoculate 200 mL LB media in a 250-mL flask supplemented with 50  $\mu$ g/mL kanamycin. Incubate the starter culture at 37°C with shaking at 250 rpm for ~14 h. Place four 6-L flasks containing 4 L M9 minimal media in the shaker overnight at 37°C to warm the culture media.
3. Prior to inoculation of the large flasks, add 50  $\mu$ g/mL kanamycin and 400  $\mu$ L 10,000 $\times$  micronutrient mix (3 nM (NH<sub>4</sub>)<sub>2</sub>Mo<sub>7</sub>O<sub>24</sub>, 30 nM CoCl<sub>2</sub>·6H<sub>2</sub>O, 400 nM H<sub>3</sub>BO<sub>3</sub>, 10 nM CuSO<sub>4</sub>·5H<sub>2</sub>O, 80 nM MnCl<sub>2</sub>·4H<sub>2</sub>O, and 10 nM ZnSO<sub>4</sub>·7H<sub>2</sub>O).
4. Use 20 mL of the overnight starter culture to inoculate each of the four large flasks.
5. At an OD<sub>600</sub> = 0.3, supplement cultures with 200 mL amino acid mixture (alanine (800  $\mu$ M), arginine (400  $\mu$ M), asparagine (400  $\mu$ M),

aspartic acid (400  $\mu$ M), cysteine (100  $\mu$ M), glutamic acid (600  $\mu$ M), glutamine (600  $\mu$ M), glycine (800  $\mu$ M), histidine (200  $\mu$ M), isoleucine (400  $\mu$ M), leucine (800  $\mu$ M), lysine (400  $\mu$ M), methionine (200  $\mu$ M), phenylalanine (400  $\mu$ M), proline (400  $\mu$ M), serine (10 mM), threonine (400  $\mu$ M), tryptophan (100  $\mu$ M), tyrosine (200  $\mu$ M), and valine (600  $\mu$ M). Take a 1 mL sample (uninduced sample) for SDS-PAGE analysis.

6. At an  $OD_{600}=0.6$ , add 1,10-phenanthroline to each of the flasks to a final concentration of 75  $\mu$ M and allow the flasks to shake for an additional 15 min.
7. Place the flasks in an ice-water bath for  $\sim 1$  h. Turn the temperature of the incubator to 18°C.
8. Induce overexpression of the desired gene by adding 200  $\mu$ M IPTG and allow expression to take place at 18°C with shaking at 180 rpm for  $\sim 18$  h.
9. On the following day, take a sample (post-IPTG induction) normalized to the number of cells in the uninduced sample for SDS-PAGE analysis.
10. Harvest the cells by centrifugation at  $7500 \times g$  for 12 min. Flash-freeze the resulting cell pellet ( $\sim 2$  g/L) in liquid nitrogen and store in liquid nitrogen until further use.
11. Isolate the apo protein using the exact method described for the holo protein.
12. Analyze the overproduction of the protein by SDS-PAGE analysis (12.5%) using the samples collected.

## 4.2 Synthesis of $^{34}\text{S}$ -Labeled Sulfide From $^{34}\text{S}$ -Labeled Elemental Sulfur

1. Perform all steps under a constant flow of nitrogen.
2. Cool a two-neck round-bottomed flask containing a magnetic stir bar to  $-78^\circ\text{C}$  in a dry ice/acetone bath.
3. Add 0.5 g elemental  $^{34}\text{S}$ -labeled elemental sulfur to the cooled flask and pass ammonia gas into the flask until approximately 75 mL liquid ammonia forms.
4. Carefully and slowly add a twofold molar excess of sodium metal to the reaction and allow it to stir for 30 min.
5. Remove the reaction from the dry ice/acetone bath and allow it to warm to room temperature.
6. Evaporate the ammonia under a stream of nitrogen gas to yield a white precipitate.
7. Add isopropanol until the blue color dissipates, indicating that the excess sodium metal is quenched.

8. Dissolve the white precipitate in 15 mL H<sub>2</sub>O.
9. Remove the solvent using a rotary evaporator to yield a solid containing both NaOH and <sup>34</sup>S-labeled sodium sulfide (Na<sub>2</sub><sup>34</sup>S).
10. Dissolve the solid in a minimum amount of water and quantify the sulfide concentration using the colorimetric sulfide analysis protocol detailed above.

### 4.3 Reconstitution of Apo NfuA Using Na<sub>2</sub><sup>34</sup>S

1. Perform all reconstitution steps under strictly anaerobic conditions.
2. Dilute Apo NfuA to 100 μM in reconstitution buffer (50 mM HEPES, pH 7.5, 300 mM KCl, 5 mM DTT, and 20% glycerol) in a total volume of 20 mL.
3. Add DTT (5 mM final concentration) in three increments every 20 min. Incubate the solution on ice for 1 h.
4. Slowly, add FeCl<sub>3</sub> (2 mM final concentration) every 5 min a total of five times. Incubate the solution for an additional 30 min.
5. Add the Na<sub>2</sub><sup>34</sup>S (2 mM final concentration) in five increments every 30 min. Check the pH of the Na<sub>2</sub><sup>34</sup>S solution before adding and neutralize the NaOH left from the synthesis procedure if necessary. Incubate the reaction on ice overnight.
6. The following day, remove protein precipitation by centrifuging the reconstituted protein at 14,000 × *g* for 15 min.
7. Concentrate the solution to ~2 mL.
8. Remove aggregates from the reconstitution procedure by size-exclusion chromatography using a HiPrep 16/60 Sephacryl HR S-200 column (GE Health Sciences) equilibrated in storage buffer (flow rate of 0.5 mL/min) and connected to an ÄKTA protein liquid chromatography system (GE Health Sciences) housed in an anaerobic chamber.
9. Collect the fractions containing the reconstituted protein and concentrate to ~2 mL. Flash-freeze in liquid nitrogen until further use. Redetermine the protein concentration using the Bradford method and perform an iron and sulfide analyses on the reconstituted protein.




## 5. STRATEGIES FOR ANALYZING PROTEIN COMPLEX FORMATION

Various methods exist to measure protein–protein interactions; however, it is often difficult to characterize interactions that are transient in nature. A nonexhaustive list of methods includes coimmunoprecipitation,

isothermal titration calorimetry (ITC), surface plasmon resonance (SPR), and analytical molecular-sieve chromatography. The benefit of using coimmunoprecipitation is the ability to pull down protein complexes *in vivo* that are biologically relevant. However, while it is possible to pull protein complexes out of the cell, the use of a cross-linking reagent is often needed to capture weak interactors. Both ITC and SPR are enticing options because they can quantitatively determine the  $K_d$  associated with the binding event in protein–protein interactions or with protein–small molecule interactions. One downfall associated with ITC is the need for large amounts of material for analysis.

Qualitatively, molecular-sieve chromatography can be used to determine if *in vitro* mixtures of proteins coelute as a complex rather than as individual components. A standard curve can be generated using a suite of standards in order to estimate the molecular weight of the complexes (Fig. 5). This method is inferior to ITC and SPR in the sense that it only provides a qualitative assessment of whether the proteins form a complex. However, it takes very small amounts of material and can be performed in a single day. Similar to coimmunoprecipitation, this method does not capture transient interactions. Following its completion, other methods can be used for a more quantitative assessment of the binding event.



**Fig. 5** Molecular sieve chromatography. (A) Chromatogram of the suite of standards used to generate the standard curve as follows: (A) Blue dextran, 2000 kDa; (B)  $\beta$ -Amylose, 200 kDa; (C) Alcohol dehydrogenase, 150 kDa; (D) Bovine serum albumin, 66 kDa; (E) Carbonic anhydrase, 29 kDa; (F) Cytochrome C, 12.4 kDa. (B) Standard curve generated from the elution volume of each of the standards corrected for the void volume of the column vs the log molecular weight of each of the standards.

## 5.1 Protocol for Performing Analytical Molecular-Sieve Chromatography

1. Equilibrate a HiPrep 16/60 Sephacryl HR S-200 column (GE Healthcare) connected to an ÄKTA protein liquid chromatography system (GE Health Sciences) housed in an anaerobic chamber with buffer (50 mM HEPES, pH 7.5, 300 mM KCl, 5 mM DTT, and 15% glycerol).
2. It is essential to ensure that the glove box is fully devoid of oxygen contamination ( $<1$  ppm  $O_2$ ) and to add fresh reductant (5 mM DTT) to buffer directly before use to prevent disulfide formation between proteins.
3. Prepare a 500  $\mu$ L sample in buffer (50 mM HEPES, pH 7.5, 300 mM KCl, 5 mM DTT, and 15% glycerol) containing the following proteins as standards: cytochrome *C* (0.3 mg/mL), carbonic anhydrase (0.3 mg/mL), bovine serum albumin (0.80 mg/mL), alcohol dehydrogenase (0.65 mg/mL),  $\beta$ -amylase (0.75 mg/mL), and Blue Dextran (0.40 mg/mL).
4. Load the mixture containing the standards onto the column and flow at a rate of 0.5 mL/min until all of the standards elute. Record the elution volume of the Blue Dextran, which indicates the experimentally determined void volume ( $V_0$ ), and the elution volumes ( $V_e$ ) of each of the additional standards, which will be used to generate a standard curve.
5. Dilute the LipA protein to make a 500  $\mu$ L sample at a concentration of 100  $\mu$ M. Load the sample onto the column and elute at a rate of 0.5 mL/min. Record its elution volume.
6. Dilute the NfuA protein to make a 500  $\mu$ L sample with a concentration of 100  $\mu$ M. Load the sample onto the column and elute at a rate of 0.5 mL/min. Record its elution volume.
7. Make a 500  $\mu$ L sample mixture containing 100  $\mu$ M LipA protein and 100  $\mu$ M NfuA protein. Incubate at room temperature for 5 min and then load onto the column and elute at a rate of 0.5 mL/min. Collect protein containing fractions and record the elution volumes of each peak observed in the chromatogram.
8. Concentrate the fractions corresponding to each protein peak observed in the chromatogram. Analyze by SDS-PAGE (12.5%) to confirm the identity of the proteins in each fraction.
9. Generate a standard curve of the Log(molecular weight) vs  $V_e/V_0$  for each of the standards.
10. Divide the experimentally determined elution volume of each of the peaks in the chromatograms of each sample by the void volume of the column ( $V_e/V_0$ ).




11. Use the linear equation from the standard curve to solve for the Log(molecular weight) of the individual proteins as well as the protein complexes.



## 6. MONITORING IRON–SULFUR CLUSTER TRANSFER USING QUANTITATIVE ANALYSIS

A strategy for accurately monitoring lipoyl cofactor production *in vitro* was not established until the early 2000s (Miller et al., 2000). This strategy was based on an indirect assay in which activation of apo-pyruvate dehydrogenase by lipoylation was coupled to the reduction of  $\text{NAD}^+$ , which could be monitored spectrophotometrically. Through this method, the rate of NADH formation was directly dependent on the amount of lipoylated pyruvate dehydrogenase. Since that time, more direct methods have been developed to quantitatively monitor the formation of the mono-thiolated intermediate and lipoyl product using liquid chromatography-mass spectrometry (LC-MS), which is described below. This method allows the activity of LipA to be quantified and the effect of NfuA on LipA's activity to be assessed (Fig. 6). Additionally, by using NfuA with isotopically labeled clusters, as detailed above, the transfer of an iron–sulfur cluster can be directly monitored. Further, including a competitive chelator such as citrate in the reaction can serve to ensure that the cluster transfer from



**Fig. 6** Time-dependent formation of the monothiolated intermediate and lipoyl product in the absence (15  $\mu\text{M}$  LipA, shown in black) and presence (400  $\mu\text{M}$  NfuA, shown in red) of NfuA.

NfuA to LipA is direct, as opposed to a diffusive mechanism wherein LipA obtains a cluster from solution (Ding, Clark, & Ding, 2004).

## 6.1 Sample Preparation for Analysis of Reaction Products by LC-MS to Assess the Effect of NfuA on LipA Activity

1. Reactions are performed anaerobically. Reactions are typically prepared in a total volume of 200  $\mu\text{L}$  and contain the following: buffer (50 mM HEPES, pH 7.5, 300 mM KCl, 10% glycerol), 600  $\mu\text{M}$  octanoyl peptide substrate (Glu-Ser-Val-[ $\text{N}^6$ -octanoyl]Lys-Ala-Ala-Ser-Asp), 25  $\mu\text{M}$  LipA, 1  $\mu\text{M}$  *S*-adenosylhomocysteine (SAH) nucleosidase (to prevent product inhibition), and 2 mM SAM. An additional and identical reaction, but containing 400  $\mu\text{M}$  NfuA, is also prepared. Importantly, DTT has been shown to be capable of nonphysiological iron-sulfur cluster transfer and should not be used in assay buffers (Vranish et al., 2015).
2. A 500  $\mu\text{L}$  quench solution containing 100 mM  $\text{H}_2\text{SO}_4$  and 20  $\mu\text{M}$  AtsA (peptide internal standard, Pro-Met-Ser-Ala-Pro-Ala-Arg-Ser-Met) is prepared, and 15  $\mu\text{L}$  is aliquoted into 11 Eppendorf tubes.
3. Initiate the reaction by adding the low-potential reductant. In LipA/NfuA assays, 2 mM dithionite is typically used. The initiation of the two reactions (with and without NfuA) can be staggered by 30 s in order to perform them simultaneously.
4. At the specified time points ( $t=0, 1, 2, 4, 8, 15, 30, 60, 90, 120$ , and 150 min), add 15  $\mu\text{L}$  of the reaction into 15  $\mu\text{L}$  of the quench solution.
5. After the reaction is complete, remove the reaction tubes from the anaerobic chamber and centrifuge for 30 min at  $14,000 \times g$ .
6. Avoiding any precipitation, pipette 25  $\mu\text{L}$  of each time point into an appropriate autosampler vial for LC-MS analysis.

## 6.2 Sample Preparation for Analysis of Reaction Products in the Presence of a Competitive Chelator or Exogenous $^{34}\text{S}$ -Labeled Sulfide to Ensure Direct Cluster Transfer From NfuA to LipA

1. Prior to performing activity assays with LipA and NfuA, incubate each protein individually with a fivefold excess sodium citrate. Assess whether the citrate has a significant effect on the Fe—S clusters of LipA and NfuA by obtaining a UV-visible spectrum of the protein before and after a 1.5 h incubation with citrate.
2. Perform a LipA only control and a LipA/NfuA activity assay exactly as described previously.

3. Set up two additional reactions containing both LipA and NfuA. In one reaction, add a fivefold excess of citrate. In the other reaction, add  $\text{Na}_2^{34}\text{S}$  to a final concentration of 1 mM.
4. When analyzing data from the LipA/NfuA reaction with citrate, assess whether the addition of this chelator inhibits product formation by LipA. Citrate can chelate metals from solution and will inhibit multiple turnover of LipA if LipA's auxiliary cluster is regenerated chemically from iron and sulfide in solution after each turnover. The absence of an effect is consistent with a direct cluster transfer from NfuA through protein–protein interactions.
5. As a second method (or if citrate chelates the cluster out of the protein of interest), analyze whether  $^{34}\text{S}$  from exogenous  $\text{Na}_2^{34}\text{S}$  included in the reaction mixture is observed in the lipoyl product. The observation of product that contains sulfur at natural abundance suggests that the attached sulfur atom did not come from solution.

### 6.3 Protocol for the Quantification of Reaction Products by LC-MS

1. Prepare 50  $\mu\text{L}$  of a 500  $\mu\text{M}$  standard in deoxygenated water that contains buffer (50 mM HEPES, pH 7.5, 300 mM KCl, and 10% glycerol), 50 mM  $\text{H}_2\text{SO}_4$ , 2 mM TCEP, 20  $\mu\text{M}$  AtsA internal standard, 500  $\mu\text{M}$  octanoyl peptide, 500  $\mu\text{M}$  6-thio-octanoyl peptide, and 500  $\mu\text{M}$  lipoyl peptide. The peptides used for quantification (octanoyl, 6-thio-octanoyl, lipoyl, AtsA) are >98% pure and can be purchased from ProImmune Custom Peptides.
2. Prepare 300  $\mu\text{L}$  of a dilution mix containing the following: buffer (50 mM HEPES, pH 7.5, 300 mM KCl, and 10% glycerol), 50 mM  $\text{H}_2\text{SO}_4$ , 2 mM TCEP, and 20  $\mu\text{M}$  AtsA internal standard.
3. Aliquot 25  $\mu\text{L}$  of the dilution mix into 11 Eppendorf tubes.
4. Perform a 1:2 serial dilution of the 500  $\mu\text{M}$  standard by mixing 25  $\mu\text{L}$  of the standard with 25  $\mu\text{L}$  of the dilution mix.
5. Continue serial diluting the standards by a factor of 2 to generate standards for quantification (500, 250, 125, 62.5, 31, 16, 8, 4, 2, 1, 0.5, and 0.25  $\mu\text{M}$ ).
6. Analyze the standards and reaction mixtures by LC-MS (Agilent Technologies 1200 system coupled to an Agilent Technologies 6410 QQQ mass spectrometer) using an Agilent Technologies Zorbax Extend-C18 column Rapid Resolution HT (4.6 mm  $\times$  50 mm, 1.8  $\mu\text{m}$  particle size) equilibrated in 98% Solvent A (0.1% formic acid, pH 2.6) and 2% Solvent B (100% acetonitrile). Detect the compounds using

electrospray ionization in positive mode ( $\text{ESI}^+$ ) with the following parameters: 340°C nitrogen gas temperature, 9.0 L/min flow rate, 40 PSI nebulizer pressure, and 4000 V capillary pressure. Apply a gradient of 2%–23% Solvent B from 0.8 to 3.5 min and maintain at 23% Solvent B until 8 min. Return to 2% Solvent B from 8 to 10 min and then allow the column to reequilibrate for 2 min before the next injection. Maintain a flow rate of 0.4 mL/min throughout the entire method.

7. Use the associated MassHunter program to quantify the decay of the octanoyl peptide substrate and formation of the 6-thio-octanoyl peptide intermediate and lipoyl product (Tables 1 and 2).

**Table 1** Parameters for LC-MS/MS Analysis of Unlabeled Reaction Intermediate and Product

Compound	Parent Ion	Fragmentator Voltage	Product Ion 1 <sup>a</sup>	Product Ion 2 <sup>a</sup>
Octanoyl peptide	932.5	240	552.3 (37)	210.1 (51)
6-Thio-octanoyl peptide	964.5	240	584.3 (37)	242.1 (51)
Lipoyl peptide	996.5	240	616.4 (37)	274.1 (51)
AtsA peptide (IS)	474.4	135	229.1 (12)	153.1 (26)
Tryptophan (IS)	188.0	130	146.1 (10)	118.0 (21)

<sup>a</sup>Respective collision energies are in parentheses (McCarthy & Booker, 2017).

**Table 2** Parameters for LC-MS/MS Analysis of <sup>34</sup>S-Labeled Reaction Intermediate and Product

Compound	Parent Ion	Fragmentator Voltage	Product Ion 1 <sup>a</sup>	Product Ion 2 <sup>a</sup>
Octanoyl peptide	932.5	240	552.3 (37)	210.1 (51)
<sup>32</sup> S-6-Thio-octanoyl peptide	964.5	240	584.3 (37)	242.1 (51)
<sup>34</sup> S-6-Thio-octanoyl peptide	966.5	240	586.3 (37)	244.1 (51)
<sup>32</sup> S/ <sup>32</sup> S Lipoyl peptide	996.5	240	616.4 (37)	274.1 (51)
<sup>32</sup> S/ <sup>34</sup> S Lipoyl peptide	998.5	240	618.4 (37)	276.1 (51)
<sup>34</sup> S/ <sup>34</sup> S-Lipoyl peptide	1000.5	240	620.4 (37)	278.1 (51)
AtsA peptide (IS)	474.4	135	229.1 (12)	153.1 (26)
Tryptophan (IS)	188	130	146.1 (10)	118.0 (21)

<sup>a</sup>Respective collision energies are in parentheses (McCarthy & Booker, 2017).



## 7. CONCLUSIONS

LipA catalyzes the attachment of sulfur atoms to an unactivated octanoyl chain in the biosynthesis of the lipoyl cofactor, an essential compound used in various multienzyme complexes in Nature and implicated in human health and disease. Refined techniques for the overproduction and isolation of Fe—S containing proteins allowed for the isolation of homogenous LipA with high activity for robust mechanistic study, although the enzyme only catalyzes a single turnover due its inactivation upon the destruction of its auxiliary cluster. Recent work has identified a pathway for reforming the auxiliary cluster to allow for multiple turnovers. Biochemical methods to study the molecular details of cofactor repair and/or reinsertion in proteins that sacrifice their Fe—S clusters will help fill outstanding gaps in our understanding of the mechanism of protein maturation.

## ACKNOWLEDGMENTS

This work was supported by the National Institutes of Health Grant GM122595 (S.J.B.) and NSF Grant MCB-1158486 (S.J.B.). S.J.B. is an investigator of the Howard Hughes Medical Institute.

## REFERENCES

- Angelini, S., Gerez, C., Ollagnier-de Choudens, S., Sanakis, Y., Fontecave, M., Barras, F., et al. (2008). NfuA, a new factor required for maturing Fe/S proteins in *Escherichia coli* under oxidative stress and iron starvation conditions. *The Journal of Biological Chemistry*, 283(20), 14084–14091.
- Beinert, H. (1978). Micro methods for the quantitative determination of iron and copper in biological material. *Methods in Enzymology*, 54, 435–445.
- Beinert, H. (1983). Semi-micro methods for analysis of labile sulfide and of labile sulfide plus sulfane sulfur in unusually stable iron-sulfur proteins. *Analytical Biochemistry*, 131(2), 373–378.
- Blaszczyk, A. J., Wang, R. X., & Booker, S. J. (2017). TsrM as a model for purifying and characterizing cobalamin-dependent radical S-adenosylmethionine methylases. *Methods in Enzymology*, 595, 303–329.
- Borcea, V., Nourooz-Zadeh, J., Wolff, S. P., Klevesath, M., Hofmann, M., Urlich, H., et al. (1999). Alpha-lipoic acid decreases oxidative stress even in diabetic patients with poor glycemic control and albuminuria. *Free Radical Biology & Medicine*, 26, 1495–1500.
- Boutigny, S., Saini, A., Baidoo, E. E., Yeung, N., Keasling, J. D., & Butland, G. (2013). Physical and functional interactions of a monothiol glutaredoxin and an iron sulfur cluster carrier protein with the sulfur-donating radical S-adenosyl-L-methionine enzyme MiaB. *The Journal of Biological Chemistry*, 288(20), 14200–14211.
- Bradford, M. M. (1976). A rapid and sensitive method for the quantitation of microgram quantities of protein utilizing the principle of protein-dye binding. *Analytical Biochemistry*, 72, 248–254.

- Cicchillo, R. M., & Booker, S. J. (2005). Mechanistic investigations of lipoic acid biosynthesis in *Escherichia coli*: Both sulfur atoms in lipoic acid are contributed by the same lipoyl synthase polypeptide. *Journal of the American Chemical Society*, 127, 2860–2861.
- Cicchillo, R. M., Iwig, D. F., Jones, A. D., Nesbitt, N. M., Baleanu-Gogonea, C., Souder, M. G., et al. (2004). Lipoyl synthase requires two equivalents of S-adenosyl-L-methionine to synthesize one equivalent of lipoic acid. *Biochemistry*, 43, 6378–6386.
- Cicchillo, R. M., Lee, K. H., Baleanu-Gogonea, C., Nesbitt, N. M., Krebs, C., & Booker, S. J. (2004). *Escherichia coli* lipoyl synthase binds two distinct [4Fe–4S] clusters per polypeptide. *Biochemistry*, 43, 11770–11781.
- Cronan, J. E. (2016). Assembly of lipoic acid on its cognate enzymes: An extraordinary and essential biosynthetic pathway. *Microbiology and Molecular Biology Reviews*, 80, 429–450.
- Ding, H., Clark, R. J., & Ding, B. (2004). IscA mediates iron delivery for assembly of iron-sulfur clusters in IscU under the limited accessible free iron conditions. *The Journal of Biological Chemistry*, 279, 37499–37504.
- Douglas, P., Kriek, M., Bryant, P., & Roach, P. L. (2006). Lipoyl synthase inserts sulfur atoms into an octanoyl substrate in a stepwise manner. *Angewandte Chemie*, 118, 5321–5323.
- Harmer, J. E., Hiscox, M. J., Dinis, P. C., Fox, S. J., Iliopoulos, A., Hussey, J. E., et al. (2014). Structures of lipoyl synthase reveal a compact active site for controlling sequential sulfur insertion reactions. *The Biochemical Journal*, 464, 123–133.
- Holmquist, L., Stuchbury, G., Berbaum, K., Muscat, S., Young, S., Hager, K., et al. (2007). Lipoic acid as a novel treatment for Alzheimer's disease and related dementias. *Pharmacology & Therapeutics*, 113(1), 154–164.
- Landgraf, B. J., McCarthy, E. L., & Booker, S. J. (2016). Radical S-adenosylmethionine enzymes in human health and disease. *Annual Review of Biochemistry*, 85(1), 485–514.
- Lanz, N. D., & Booker, S. J. (2015). Auxiliary iron-sulfur cofactors in radical SAM enzymes. *Biochimica et Biophysica Acta*, 1853(6), 1316–1334.
- Lanz, N. D., Grove, T. L., Gogonea, C. B., Lee, K. H., Krebs, C., & Booker, S. J. (2012). RlmN and AtsB as models for the overproduction and characterization of radical SAM proteins. *Methods in Enzymology*, 516, 125–152.
- Lanz, N. D., Pandelia, M. E., Kakar, E. S., Lee, K. H., Krebs, C., & Booker, S. J. (2014). Evidence for a catalytically and kinetically competent enzyme-substrate cross-linked intermediate in catalysis by lipoyl synthase. *Biochemistry*, 53, 4557–4572.
- McCarthy, E. L., & Booker, S. J. (2017). Destruction and reformation of an iron-sulfur cluster during catalysis by lipoyl synthase. *Science*, 358(6361), 373–377.
- McLaughlin, M. I., Lanz, N. D., Goldman, P. J., Lee, K. H., Booker, S. J., & Drennan, C. L. (2016). Crystallographic snapshots of sulfur insertion by lipoyl synthase. *Proceedings of the National Academy of Sciences of the United States of America*, 113, 9446–9450.
- Miller, J. R., Busby, R. W., Jordan, S. W., Cheek, J., Henshaw, T. F., Ashley, G. W., et al. (2000). *Escherichia coli* LipA is a lipoyl synthase: In vitro biosynthesis of lipoylated pyruvate dehydrogenase complex from octanoyl-acyl carrier protein. *Biochemistry*, 39, 15166–15178.
- Pandelia, M. E., Lanz, N. D., Booker, S. J., & Krebs, C. (2015). Mossbauer spectroscopy of Fe/S proteins. *Biochimica et Biophysica Acta*, 1853(6), 1395–1405.
- Py, B., Gerez, C., Angelini, S., Planel, R., Vinella, D., Loiseau, L., et al. (2012). Molecular organization, biochemical function, cellular role and evolution of NfuA, an atypical Fe-S carrier. *Molecular Microbiology*, 86(1), 155–171.
- Reed, L. (1974). Multienzyme complexes. *Accounts of Chemical Research*, 7, 40–46.
- Reed, L. J., Busk, B. G. D., Gunsalus, I. C., & Schnakenberg, G. H. F. (1951). Crystalline  $\alpha$ -lipoic acid: A catalytic agent associated with pyruvate dehydrogenase. *Science*, 114, 93.
- Reed, L. J., Gunsalus, I. C., Schnakenberg, G. H. F., Soper, Q. F., Boaz, H. E., Kern, S. F., et al. (1953). Isolation, characterization and structure of  $\alpha$ -lipoic acid. *Journal of the American Chemical Society*, 75, 1267–1270.

- Vranish, J. N., Russell, W. K., Yu, L. E., Cox, R. M., Russell, D. H., & Barondeau, D. P. (2015). Fluorescent probes for tracking the transfer of iron-sulfur cluster and other metal cofactors in biosynthetic reaction pathways. *Journal of the American Chemical Society*, 137(1), 390–398.
- Yi, X., & Maeda, N. (2005). Endogenous production of lipoic acid is essential for mouse development. *Molecular and Cellular Biology*, 25, 8387–8392.
- Zhao, S., Miller, J. R., Jiang, Y., Marletta, M. A., & Cronan, J. E., Jr. (2003). Assembly of the covalent linkage between lipoic acid and its cognate enzymes. *Chemistry & Biology*, 10, 1293–1302.

## FURTHER READING

- Frey, P. A., & Booker, S. J. (2001). Radical mechanisms of S-adenosylmethionine-dependent enzymes. *Advances in Protein Chemistry*, 58, 1–45.
- Snell, E. E. (1992). Microorganisms in vitamin and biofactor research. *Journal of Nutritional Science and Vitaminology*, 38, 34–39. Spec No:34–9.
- Snell, E. E. (1993). From bacterial nutrition to enzyme structure: A personal odyssey. *Annual Review of Biochemistry*, 62, 1–27.