

Cryogenic transmission electron microscopy for materials research

David W. McComb, Jeffrey Lengyel, and C. Barry Carter

David W. McComb, The Ohio State University, USA; mccomb.29@osu.edu

Jeffrey Lengyel, Thermo-Fisher, USA; jeffrey.lengyel@thermofisher.com

C. Barry Carter, Storrs, CT and CINT, Sandia National Labs, USA;
cbarrycarter@gmail.com

Cryogenic transmission electron microscopy is simply transmission electron microscopy conducted on specimens that are cooled in the microscope. The target temperature of the specimen might range from just below ambient temperature to less than 4 K. In general, as the temperature decreases, cost increases, especially below -77°C when liquid He is required. We have two reasons for wanting to cool the specimen—improving stability of the material or observing a material whose properties change at lower temperatures. Both types of study have long histories. The cause of excitement in this field today is that we have a perfect storm of research activity—electron microscopes are almost stable with minimal drift (we can correct what drift there is), we can prepare specimens from the bulk or build them up, we have spherical-aberration-corrected lenses and monochromated beams, we have direct-electron-detector cameras, and computers are becoming powerful enough to handle all the data we produce.

Keywords: transmission electron microscopy (TEM). Spectroscopy, water, polymer, nanostructure, scanning transmission electron microscopy (STEM)

Introduction

Transmission electron microscopy (TEM) today is the essential tool for materials characterization when site-specific chemical or structural detail is required.^{1,2} If you want to know the chemistry at a crack tip or the uniformity of a coating on a nanoparticle, then you must use TEM or scanning transmission electron microscopy (STEM). (S)TEM can be used to characterize all materials, from the hardest ceramics and metals to the softest polymers and the smallest nanoparticles.

It can also be applied to living material, almost invariably when it is no longer living. By “living,” we mean human, animal, and vegetable, and by “no longer” we could mean after recent surgery or after mummification³ or fossilization.⁴

The TEM research community has tended to separate into three groups—the physical scientists, the biologists, and the instrumentalists. Modeling and software specialists tend to segregate into these same groups. Recent developments in TEM and STEM have begun to see a convergence of these different fields; the use of cryogenic techniques is one aspect where these groups have more in common than differences. The commonalities can be the instruments, the many ways we now have to prepare specimens for the TEM, or even the ways that we acquire, manipulate, and analyze experimental data.

Insert Figure 1 Comparison of the macroscopic appearance of cryotemplated and PGX CNC aerogels in dry and wet states. The difference is clear in the macroscopic image but cryo-TEM will show the full picture.⁵

Method of the year and a Nobel Prize

In 2017, the Nobel Prize for chemistry was awarded to Jacques Dubochet, Joachim Frank and Richard Henderson “for developing cryo-electron microscopy for the high-resolution structure determination of biomolecules in solution.”⁶ In this application of cryogenic (cryo)-TEM, wherein proteins are imaged, the temperature of the sample must be maintained at liquid-nitrogen temperatures or below. While the low temperature not only confers radiation damage protection to samples, it also is necessary to maintain the sample in the vitreous state where proteins adopt near-native conformations. Modern automated cryo-TEMs are designed specifically to operate at liquid-nitrogen temperature.

In awarding the Nobel Prize, the Royal Swedish Academy of Sciences explained that cryo-EM has “taken biochemistry into a new era.”⁷ (Microscopists were initially so excited to get some recognition that they forgot about the missing T, but we will include it.) The journal *Nature Methods* had actually named “single-particle cryo-electron microscopy (cryo-EM)” the “Method of the Year for 2015,” stating that it meant the end of “blob-ology.”⁸ The announcement was accompanied by a historical perspective,⁹ which recognized the pioneering work

by Taylor and Glaeser,¹⁰ and by a commentary asking “how good can cryo-EM become?”¹¹

Today’s success of cryo-TEM builds on the independent work of the three Nobel Laureates, but three much more recent developments were essential before its potential could be realized. We now have a new class of direct electron detector cameras, much more stable specimen stages, and for the highest resolution, an aberration-corrected objective lens. The direct electron detector camera makes every electron count by counting every electron. Very stable stages are now available because of the attention being paid to mechanical and thermal stability and such small amounts of drift can be corrected on the fly. Aberration-corrected microscopes are now available from several manufacturers thanks to developments by Rose, Haider, Urban, and Krivanek.^{12,13} While aberration correction has found limited application in cryo-TEM for structural biology, at least to date, the ultrastable electron optics necessary to support aberration correction are a key contributor to the success of cryo-TEM. Cryo-TEM is now the method of choice for characterization in structural biology, in many cases replacing X-ray diffraction, and the plethora of newly installed instruments worldwide is generating data faster than can be fully analyzed.

The rest of materials

We were then led to ask “what is the impact of cryo-TEM on materials research where the materials are not biological?” This question is the subject of this collected set of articles in this issue of *MRS Bulletin*. We lead off with an article¹⁴ discussing the instrumentation since it is the development of instruments that has made this revolution possible for biochemists. The advantages and disadvantages of TEM and STEM for biological and physical sciences research are reviewed and detector technology, both current and emerging, are discussed.

We examine the two types of material that are most likely to be impacted by cryo-TEM. These are, of course, materials that have much in common with the biological materials, namely soft materials¹⁵ and very soft materials¹⁶. The latter, liquids, are much too “soft,” so we need to use cryo-techniques to make them stay where we put them in the microscope. In both cases, it is hoped that freezing the

specimen will also stop, or at least decrease, changes that might occur over time, especially when the specimen is subjected to the electron beam. But changes will inevitably occur as they do in all TEM studies, so we include a critical consideration of what might happen when the electron beam passes through the specimen.¹⁷

Another article¹⁸ addresses another of the strengths of TEM; the TEM is an analytical machine that not only produces images, but also provides diffraction data¹⁰ and information on the chemistry, both composition and bonding. A contribution by Minor et al.¹⁹ considers the impact that cryo-EM could have in the field of quantum materials where investigation and understanding of low-temperature phenomena is critical to future developments.

It might be anticipated that the greatest impact for cryo-TEM in materials research would be in the field of polymer science or possibly in the study of composite materials that include polymers and materials that, like biological samples, are wet in real life, or may include a wet component.

Making specimens of composite materials has always been a challenge, which is why bottom-up specimen preparation is always preferred over top-down; the interest in nanoparticles has been a boon for TEM because we already have the specimen thin enough for electron transmission.

Insert Figure 2 Cryo-SEM image of the gel network formed by PVA cross-linked by borax. This hydrogel is used in a soft robotics application as described by Preller et al.²⁰

Liquid films on small objects,²¹ or even in thin layers,²² are technologically extremely important. But if you have a “bulk” sample and need to thin it, then you need to know where to go to make that specimen. Once you have the specimen, which microscope should you use? The first choice (the one you have immediate access to) may not be the best choice. Fortunately, user facilities are developing across the US and internationally, where one can collaborate to use the optimum TEM once you know what that is. The most popular collaborator for any TEM expert is a researcher who has already prepared the specimen, and the most popular collaborator for the researcher with the specimen is the TEM expert who will to examine it.

The article by Williams et al.¹⁴ describes some of the instruments available and considers how best to utilize them for different types of experiment. Not all studies require an aberration-corrected instrument. Sometimes 100 kV will be better than 300 kV depending on the material and possibly the detector. However, the flexibility to use a range of detectors and methods often comes at the cost of instrument stability.

If your radiation-sensitive sample forms nanocrystals

It is now possible to utilize a new methodology called micro-electron diffraction (MicroED) to determine the three-dimensional structure of radiation-sensitive samples such as proteins, small molecules, or nanomaterials from nanocrystals. While diffraction methods in TEM are common when applied to nanomaterials, this is a new variation of cryo-TEM.²³ In MicroED, the cryo-TEM is operated in diffraction mode, and crystals are often on the order of 1 μm in size. For protein structure determination, crystals are often a billionth of the volume typically used for x-ray crystallography. The TEM used is essentially the same as traditional cryo-TEMs, aside from a few minor modifications, and the addition of a specialized complementary metal oxide semiconductor (CMOS)-based detector optimized for continuous acquisition of extremely low-dose diffraction-pattern collection. Continuous-rotation diffraction datasets are processed using existing X-ray crystallographic software packages. Given the minor modification needed to enable MicroED, many existing cryo-TEMs can be retrofitted to perform this technique, and it is expected that there will be dramatic growth in this research methodology in the very near future. The article by Nannenga and Gonen²¹ in this issue discusses MicroED.

If it's not bio

The issue includes two articles on the materials most obviously benefitting from cryo-TEM. We consider cryo-TEM of soft materials¹⁵ and materials that would contain liquids¹⁶ if they were at room temperature. In both cases, the materials have much in common with biomaterials—the aim of cryo-temperatures is to stop, or at least minimize, any movement or changes in the specimen, so the cooling

(freezing) needs to occur as quickly as possible. When water is present, cooling has to be so fast that the water cannot crystallize—it must remain in the vitreous state. We also want no expansion to occur since this might distort the specimen. Similarly, the formation of ice crystals would add confusion to the interpretation in addition to changing the specimen.

Insert Figure 3 A cryo-TEM image of water soluble magnetite.²⁴

Clearly then, we do not want the specimen to change before we look at it in the TEM, but we also do not want it to change when we examine it. This brings us to consider how the electron beam might change the specimen, and what voltage to use. Basically, the electron beam will change the specimen; the uncertainty is only how fast it will change or, in other words, how long we have to make observations. It should be remembered that high-voltage TEMs were in fact extensively used in the 1960s to simulate damage occurring in nuclear reactors. The changes can be caused either by direct impact of electrons (energy and momentum transfer) or by interaction between charged particles—the electrons in the beam and atomic nuclei in the specimen. In the first case, it is worse if the accelerating voltage is high, while in the latter, if the voltage is low. These differences mean that the choice of TEM may not be immediately obvious for a particular specimen. Only recently has quantitative analysis of beam damage started to become a reality for organic films.²⁵ The article in this issue on solid–liquid interfaces has many implications—the liquid is not necessarily water. So the technique has enormous application to energy materials such as batteries, and fields such as corrosion.²⁶

Our penultimate topic emphasizes why the TEM is unique and aims to encourage users to become more familiar with the other capabilities of the instrument. Even if you know that your specimen only contains carbon, oxygen, and hydrogen, the TEM can give you information on the bonding; it may also tell you about contamination that you did not know was present. The TEM gives two principal approaches to chemical analysis, either x-ray energy dispersive spectroscopy (XEDS) or electron energy-loss spectroscopy (EELS); if you are making measurements, you would use spectrometry instead. Traditionally, XEDS

has favored heavy elements and EELS has favored light elements, but improved detectors and software have begun to blur these distinctions. Similarly, diffraction used to be just for TEM not STEM, but that too is changing, again partly because of new detectors (in this case the direct electron detector cameras) and partly because of improving software.²⁷ When we say software, we are not just thinking about interpreting the spectrum or diffraction pattern; today, we may be collecting so much data that no person will ever have the time to look at it all. This raises other questions such as does your institution really have the capability to store data as you claimed in your data plan and would anyone else ever have time to access it if they did. Here, the non-bio materials researchers have exactly the same challenges as the biological and biochemical cryo-TEM users.

There has long been great interest in using cryo-TEM to study materials that have different properties at cryogenic temperatures. For example, the late A. Tonomura and his colleagues performed classic experiments verifying the Aharonov–Bohm effect, observing the movement of magnetic flux lines (flux vortices) in niobium at temperatures as low as 4.2 K, and especially watching changes occur as the specimen warmed from 4.5 K and 15 K.²⁸ Such cryogenic work requires a stage that is cooled by liquid helium. These experiments are difficult because of the temperature and the need for a Lorentz lens, but the information is not achievable by any other means.²⁹ Unfortunately, TEM using liquid helium quickly becomes expensive. What physics might be revealed by an ultrastable cryo-TEM operating at liquid helium temperatures.¹⁹ Arguably, the greatest impact would be in driving quantum materials research that has become the focus of physics and materials communities across the world.

It won't cool off

We can draw a few conclusions regarding cryo-TEM and what it can do for non-bio materials science. The microscopes are essentially the same—for both groups, the TEM is an electron-optics column with an electron source and everything controlled by computers. No one ever “sees” the sample—we all form images of the specimen and then interpret the images. Both groups are interested in obtaining information about all three dimensions (not just a two-dimensional

projection) and will want to know the chemistry of the specimen. Using new direct electron-detector cameras, both groups are faced with storing, processing, and retrieving more data in a week than what previous generations of TEM users produced in a lifetime. In science and life, most things happen at interfaces. Most materials we use are composite materials, ranging from electronic devices to reinforced concrete. For medical implants, the contact to the body is critical. Targeted drug delivery uses carrier particles to transport the drug to the target.³⁰ These examples emphasize the developing links between bio- and non-bio materials. So it is not just what can we learn from the success of bio-cryo-TEM, but how can we improve the links between different fields of cryo-TEM.

Acknowledgments

C.B.C. acknowledges support from the National Science Foundation Award Nos. 1820565 and DMR-1710468. He conducts his research in CINT, which is an Office of Science User Facility operated for the US Department of Energy (DOE). Sandia National Labs is managed and operated by National Technology and Engineering Solutions of Sandia, LLC., a wholly owned subsidiary of Honeywell International, Inc., for the US DOE's NNSA contract DE-NA-0003525.

References

1. D.B. Williams, C.B. Carter, *Transmission Electron Microscopy: A Textbook for Materials Science*, 2nd ed. (Springer, New York, 2009).
2. C.B. Carter, D.B. Williams, *Transmission Electron Microscopy: Diffraction, Imaging, and Spectrometry* (Springer, Heidelberg, 2016).
3. M.A. Pabst, F. Hofer, *Am. J. Phys. Anthropol.* **107**, 1 (1998).
4. E.J. Javaux, A.H. Knoll, M.R. Walter, *Geobiology* **2**, 121 (2004).
5. D.A. Osorio, B. Seifried, P. Moquin, K. Grandfield, E.D. Cranston, *J Mater Sci* **53**, 9842 (2018)
6. The Royal Swedish Academy of Sciences, "The Nobel Prize in Chemistry 2017," <https://www.nobelprize.org/prizes/chemistry/2017/press-release/>. Press release.
7. Royal Swedish Academy of Sciences, "They captured life in atomic detail" (2017).

8. Editor, *Nat. Methods* **13**, 1 (2015).
9. E. Nogales, *Nat. Methods* **13**, 24 (2015).
10. K.A. Taylor, R.M. Glaeser, *Science* **186**, 1036 (1974).
11. R.M. Glaeser, *Nat. Methods* **13**, 28 (2015).
12. M. Haider, S. Uhlemann, E. Schwan, H. Rose, B. Kabius, K. Urban, *Nature* **392**, 768 (1998).
13. O.L. Krivanek, N. Dellby, A.R. Lupini, *Ultramicroscopy* **78**, 1 (1999).
14. R.E.A Williams, D.W. McComb, S. Subramanian, MRS Bulletin 44(12) 2019
15. J. Watt, D.L. Huber, P.L. Stewart, MRS Bulletin 44(12) 2019
16. M.J. Zachman, N. de Jonge, R. Fischer, K. Jungjohann, D. Perea, MRS Bulletin 44(12) 2019
17. C.J. Russo, R.F. Egerton, MRS Bulletin 44(12) 2019
18. B.L. Nannenga, T. Gonen, MRS Bulletin 44(12) 2019
19. A. Minor, P. Denes, D.A. Muller, MRS Bulletin 44(12) 2019
20. Preller, G. Runge, S. Zellmer, D. Menzel, S.A. Saein, J Peters. A. Raatz, B. Tiersch, J. Koetz, G. Garnweitner, *J Mater Sci* **54**, 4444 (2019)
21. J.S. Du, P.-C. Chen, B. Meckes, E.J. Kluender, Z. Xie, V.P. Dravid, C.A. Mirkin, *J. Am. Chem. Soc.* **140**, 7213 (2018).
22. C.B. Carter, P.M. Williams, *Philos. Mag.* **26**, 393 (1972).
23. B.L. Nannenga, T. Gonen, *Nat. Methods* **16**, 369 (2019).
24. J. Watt, D.L. Huber, and the Anschutz CryoEM facility staff, unpublished
25. Z. Leijten, A.D.A. Keizer, G. de With, H. Friedrich, *J. Phys. Chem. C* **121**, 10552 (2017).
26. M.J. Zachman, Z. Tu, S. Choudhury, L.A. Archer, L.F. Kourkoutis, *Nature* **560**, 345 (2018).
27. L.A. Earl, V. Falconieri, J.L. Milne, S. Subramaniam, *Curr. Opin. Struct. Biol.* **46**, 71 (2017).
28. A. Tonomura, *Rev. Mod. Phys.* **59**, 639 (1987).

29. G. Berruto, I. Madan, Y. Murooka, G.M. Vanacore, E. Pomarico, J. Rajeswari, R. Lamb, P. Huang, A.J. Kruchkov, Y. Togawa, T. LaGrange, D. McGrouther, H.M. Ronnow, F. Carbone, *Phys. Rev. Lett.* **120**, 117201 (2018).
30. K.Y. Win, S.S. Feng, *Biomaterials* **26**, 2713 (2005).

Figure captions

(Figures are copyright other publishers—references are given)

Figure 1 Comparison of the macroscopic appearance of cryotemplated and PGX CNC aerogels in dry and wet states. The difference is clear in the macroscopic image but cryo-TEM will show the full picture.⁵

Figure 2 Cryo-SEM image of the gel network formed by PVA cross-linked by borax. This hydrogel is used in a soft robotics application as described by Preller et al.²⁰

Figure 3 A cryo-TEM image of water soluble magnetite.²⁴

David William McComb is a professor of materials science and engineering and director of The Center for Electron Microscopy and Analysis at The Ohio State University. He received his BS degree in chemistry from the University of Glasgow, Scotland, and his PhD degree in physics at the University of Cambridge, UK. His research focuses on the development and application of scanning transmission electron microscopy and electron energy-loss spectroscopy as a subnanometer-scale probe of chemistry, structure, and bonding. He is a Fellow of the Royal Society of Chemistry, the Institute of Materials, and the Microscopy Society of America. McComb can be reached by email at mccomb.29@osu.edu.

Jeff Lengyel is the Principal Scientist for Electron Microscopy Life Sciences for Thermo Fisher Scientific, formerly FEI company. His previous roles include Product Marketing Manager for Structural Biology, onsite Team Lead for the FEI/NIH Living Lab for Structural Biology, and Senior Titan Krios Applications Specialist. Jeff received his Ph.D. in Biochemistry from the Univ. of Cambridge (UK) as a Univ. of Cambridge/National Institutes of Health, Health Sciences Scholar. As per this scholarship he was co-mentored by Prof. Richard Perham (Cambridge) and Dr. Sriram Subramaniam (NIH). For his doctoral research Jeff studied the Pyruvate dehydrogenase using a combination of CryoEM and Biophysical techniques. Lengyel can be reached by email at jeffrey.lengyel@thermofisher.com.

C. Barry Carter is a CINT Distinguished Affiliate Scientist at Sandia National Labs in Albuquerque, New Mexico, a Visiting Academic at the University of Manchester in the UK and a Research Professor at UConn. He was previously on the Faculty of Cornell University (1979-1991) and the University of Minnesota (1991-2007). He has a D.Phil. from Oxford University and an Sc.D. from Cambridge University. He co-authored textbooks on Transmission Electron Microscopy and Ceramic Materials. He is a Fellow of MRS, AAAS, ACerS,

RMS, MSA and MAS. His current research focuses on phase-change materials, crystal dislocations, and reactions between Li/Na and layer materials. Carter can be reached by email at cbarrycarter@gmail.com.

Figure:

Water soluble magnetite (Fe_3O_4) nanoparticles: conventional TEM could not be used to confirm these observations due to drying artifacts, however cryo-TEM revealed that these higher order assemblies were present in solution.
(Courtesy Watt and Huber, unpublished.)